**Edwards**

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Manufacturer   
Edwards Lifesciences LLC  
One Edwards Way  
Irvine, CA 92614 USA

**EC REP**  
Edwards Lifesciences Services GmbH  
Edisonstr. 6  
85716 Unterschleißheim  
Germany

02/19  
10020191002 A  
© Copyright 2019, Edwards Lifesciences LLC  
All rights reserved.

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## DUAL STAGE VENOUS RETURN CANNULA

### Instructions for Use

**CAUTION:** Federal law (U.S.A.) restricts this device to sale by, or on the order of, a physician.

#### DESCRIPTION

Edwards Lifesciences venous return cannulae are soft cannulae with encapsulated steel wire reinforcement. Venous return cannulae are offered with various tip styles, hole patterns, French sizes, lengths, straight or angled tubing, and two sizes of connectors. Some styles are offered with obturators. Each Edwards Lifesciences device is packaged sterile and non-pyrogenic in a sealed, peel-type pouch.

#### INDICATIONS FOR USE

Dual drainage venous return cannulae are indicated for venous cannulation so that extracorporeal circulation of the venous blood to a heart-lung machine may be achieved, for a duration of  $\leq 6$  hours.

Venous cannulae may be used in pediatric populations or adult populations based on flow rate requirements and patient anatomy. Please see labeling for maximum flow rate information.

#### CONTRAINDICATIONS FOR USE

Single tube venous drainage from the right atrium should not be used in procedures where air from the procedure site may enter the cannula.

This device is not intended for use other than as indicated and should not be used when any physical impairment would contraindicate its use.

#### WARNINGS AND PRECAUTIONS

Supplied sterile and non-pyrogenic in undamaged package. Do not use if device shows signs of damage (i.e., cuts, kinks, crushed areas, leakage), or if package is damaged or open as this may indicate compromised sterility and/or product damage.

**FOR SINGLE USE ONLY.** This device is designed, intended, and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE THIS DEVICE. There are no data to support the sterility, non-pyrogenicity, and functionality of the device after reprocessing. Such action could lead to illness or an adverse event, as the device may not function as originally intended.

This device is intended for short-term use only ( $\leq 6$  hours). There are no data to support the function and performance of the device beyond six hours of use.

Products are known to contain phthalates, which may be found in device materials containing plasticizers such as DEHP and BBP. High exposure to such phthalates during medical treatments in children and pregnant or nursing women may raise a concern. A review of available data and literature supports the conclusion that the benefits outweigh the overall residual risk.

When using small diameter cannulae, do not exceed the maximum pressure limit while increasing flow. A rapid pressure rise may result.

If increased resistance is felt at any time upon insertion or removal of the cannula, stop and investigate the cause before continuing. Inability to easily advance or remove these devices may indicate vascular disease or injury. Closely examine the device position within the vessel using fluoroscopy and/or ultrasound imaging prior to proceeding.

Use the largest cannula possible, consistent with good surgical technique, to ensure adequate venous drainage.

Line pressure exceeding acceptable clinical limits may result from incorrect cannula tip position and/or restricted tip patency. Confirm proper placement to minimize risk of injury.

Avoid over-insertion of the cannula to reduce risk of cannula tip impingement against the vessel walls which may lead to decreased drainage and vessel damage.

Overinsertion of the cannula may result in proximal holes entering the inferior vena cava. This may lead to diminished venous return. Cannula markings and/or palpation is recommended to ensure proper cannula placement.

Prior to initiating cardiopulmonary bypass, completely expel air from system. Failure to eliminate air from the system could result in an air embolus.

Surgical entrance or introduction of air into the right heart may result in the partial or total loss of venous drainage to the extracorporeal circulation.

Securely tie-band the connector to cannula junction prior to initiating bypass to protect against tubing slippage.

Wire-reinforced cannulae should be clamped in the non-reinforced section located at the connector end since clamping of the reinforced section may produce permanent cannula deformation, thereby impeding flow through the cannula and risking puncture or tearing of the cannula.

Ensure proper levels of anticoagulant therapy are maintained prior to insertion of the cannula and throughout cardiopulmonary bypass, to prevent thrombus formation on or within the cannula, and in the blood stream.

Proper surgical procedures and techniques are the responsibility of the medical professional. Described procedures are provided for informational purposes only. Each physician must determine the appropriate use of this device for each patient based on medical training, experience, the type of procedure employed, and the benefits and risks associated with device use.

Dispose of used product in accordance with established hospital protocols for biohazards to minimize risk of exposure.

Graphs of pressure-flow relationships are provided at the end of these instructions.

#### DIRECTIONS FOR USE

1. Using standard surgical technique, place cannula so that drainage of the superior vena cava and right atrium occurs through proximal side holes and drainage of the inferior vena cava occurs through cannula tip hole(s).

Note: The depth markings on dual stage venous cannulae are intended to be visual rather than numerical indicators that provide reference for positioning the drainage basket within the right atrium.

2. Secure cannula. Connect to blood tubing for extracorporeal circulation.

3. Upon completion of procedure, remove cannula and repair insertion site. Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

**WARNING:** When using small diameter cannulae, do not exceed the maximum pressure limit while increasing flow. A rapid pressure rise may result.

Edwards, Edwards Lifesciences, and the stylized E logo are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

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