ASAHI

Neurovascular Guide Wire

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

ASAHI Neurovascular Guide Wire

INSTRUCTIONS FOR USE

Read these instructions carefully before use and observe the Indications for Use, Contraindications, Warnings, Precautions, Malfunctions and adverse effects and How to Use sections in the Instructions for Use. Failure to do so may result in complications, including serious injury to the patient or death.

Descriptions

This guide wire has a coil-type distal end. The coil is wholly radiopaque to facilitate selection of the blood vessel and confirmation of the position of the guide wire's distal end by fluoroscopy. The distal end can be shaped. The core shaft surface is coated with hydrophilic coating. For specifications of the respective product, refer to the product label.

Indications for Use

This guide wire is intended to be used in the neuro vasculature to facilitate the placement and exchange of therapeutic devices such as cerebral catheters during neuroradiology. **This guide wire is intended for use only in the neuro vasculature.**

Contraindications

None

Warnings

- This guide wire is for single use only. Do not reuse or resterilize. If reused or resterilized, the performance or quality of this guide wire may be compromised and there is a risk of complications, including infection.
- Do not use the guide wire after the expiration date indicated on the label. Discard any device that exceeds the expiration date.
- This guide wire must be used only by a physician who is fully trained in neuroradiology treatment.
- This guide wire must be used in an institution where emergency surgical operation can be performed immediately. [If an emergency surgical operation is unavailable, in the worst case, life-threatening events may occur.]
- Do not use in areas of vessel that are not or cannot be visualized.
- Do not modify this guide wire regardless of the reason.
- The coil section is especially fragile, so do not bend or pull it more than necessary. Otherwise, the guide wire may be damaged.
- Do not use a damaged device. Using a damaged device may result in blood vessel damage and/or inaccurate torque response. Injury to the patient may result.
- Always advance and withdraw the guide wire slowly.
- Observe movement of this guide wire in the vessels. Before this guide wire is moved or torqued, the tip movement should be examined and monitored under fluoroscopy. Do not move or torque the guide wire without observing corresponding movement of the tip. Otherwise, the guide wire may be damaged and/or vessel trauma may occur. In addition, ensure that the distal tip of this guide wire and its location in the vessel are visible during manipulations of the guide wire.
- Never push, auger, withdraw, or press this guide wire with enough force to feel resistance. Pressing or pushing this guide wire against resistance may cause damage and/or tip separation of this guide wire or direct damage to a vessel. Resistance may be felt and/or observed under fluoroscopy by noting any buckling of the tip of the guide wire. If the prolapse of the guide wire tip is observed, do not allow the tip to remain in a prolapsed position. Otherwise damage to the guide wire may occur. Determine the cause of resistance under

fluoroscopy and take any necessary remedial action.

- If any resistance is felt due to spasm or the guide wire being bent or trapped while operating the guide wire in the blood vessel or removing it, do not move or torque the guide wire. Stop the procedure. Determine the cause of resistance under fluoroscopy and take appropriate remedial action. If the guide wire is moved excessively, it may be damaged including separation or the like, which may cause blood vessel injury or result in fragments being left inside the vessel.
- When torquing this guide wire inside the blood vessel, do not torque continuously in the same direction. It may be damaged including separation or the like, which may injure the blood vessel or leave fragments inside the vessel. When torquing the guide wire, rotate it clockwise and counterclockwise alternately. Do not exceed two rotations (up to 720°) in the same direction.
- Do not push the guide wire more than necessary to advance the tip through the narrowed part of the vessel. (For example, do not push the guide wire when the distal tip of the guide wire is bent by the force of manipulation.) After crossing the targeted area, do not roughly twist, push or pull the guide wire. If the device is moved excessively, it may be damaged including separation or the like, which may injure the blood vessel or leave fragments inside the vessel.
- When inserting or removing other devices over the guide wire, flush the concomitant device with heparinized saline or other suitable solution to prevent air from entering it. Perform exchange of this guide wire carefully to prevent air entry and/or trauma. When reintroducing the guide wire, confirm that the tip is free within the vessel lumen and is not against the vessel wall. Failure to do so may result in vessel trauma when the guide wire is removed. Use the radiopaque marker of the interventional device to confirm position.
- Free movement of the guide wire within the interventional device is an important feature of a steerable guide wire system because it gives the user valuable tactile information. Test the system for any resistance prior to use. Adjust or replace the hemostatic valve with an adjustable valve if it is found to inhibit the guide wire movement.
- Before use, inspect carefully and confirm all devices and packages are undamaged.
- Before use, confirm that the guide wire is compatible with the interventional device to be used.
- Before use, confirm that the flexibility, shape, and size of the distal end of this guide wire are compatible with the procedure.
- The distal end excels in vascular selectivity. However, if it is advanced into a thin peripheral area of the vessel, the penetrating force in the vessel becomes stronger, and there will be a higher risk of vessel damage such as vessel perforation. Therefore, when manipulating the guide wire or having a micro catheter or the like follow it, take extreme care. In particular, in a procedure for collecting or suctioning blood clots for acute cerebrovascular occlusion, the guide wire manipulation becomes a blind manipulation to an occluded peripheral area for which there is no anatomical information of the vessel, therefore, more careful manipulation will be required. Do not perform catheter manipulation with the distal end of the guide wire being secured to the collateral since it will increase the risk of vessel perforation.
- Since the guide wire is provided with hydrophilic coating in order to increase the lubricity, handle it with extreme care.
- A manipulation through the lumen of the stent body using the guide wire may damage the stent or cause breakage such as separation of the guide wire. Therefore, if it is necessary to perform such a manipulation, take extreme care.
- If it is suspected that the guide wire manipulation cannot be continued, do not forcibly manipulate the guide wire, but pay a great deal of attention to any malfunctions or adverse effects, and carefully remove the whole system. If it is judged that removal of the whole system may cause malfunctions or adverse effects, immediately abandon the neuroradiology, and perform an emergency surgery procedure based on judgment of a physician.

- During the surgery, provide appropriate anticoagulation while considering the condition of the patient.
- Do not use the guide wire for pregnant patients or patients who may be pregnant (X-ray may affect the fetus).
- Do not use the guide wire for patients who cannot adapt to emergency surgery or patients who have experienced a serious allergic reaction to drugs required for the procedure such as a contrast medium.
- Do not use the guide wire for patients for whom a physician judges it inappropriate to use.
- Be sure to refer to the package insert and instructions for use of the medical equipment and drug to be used with the guide wire.

Precautions

- If the package is opened or damaged, do not use the guide wire. Do not open the package until just prior to use. Use aseptic technique in handling and using the guide wire.
- Contraindications, warnings, precautions, and intended uses of therapeutic devices that are compatible with the guide wire are described in the Instructions for Use supplied with the respective therapeutic devices. Before using the guide wire with other therapeutic devices (Sheath introducer, Shaping device, Neuroradiology guide wire, Neuroradiology guiding catheter, Micro catheter, Detachable coil, Neuroradiology dilatation catheter and Stent), read the Instructions of those devices to ensure the other devices are compatible with the guide wire. Ensure you choose a guide wire suitable for the therapy and that its use is consistent with the contraindications, warnings, precautions, and Instructions for Use of both the other devices and the guide wire.
- Guide wires are delicate instruments and should be handled carefully. When taking the guide wire out of the holder tube, do not handle the device roughly or pull it out abruptly.
- Inspect the guide wire carefully for bends, kinks, or other damage prior to use and whenever possible during the procedure.
- Never use metallic needles or metallic sheaths for insertion and withdrawal of this guide wire. Otherwise, the surface of this guide wire may be damaged significantly.
- Never use this guide wire with a metallic dilator, catheter, etc. whose metallic part may come into direct contact with the surface of this guide wire.
- When shaping the distal end, use the minimum force needed so that the coil is not damaged. Inspect the coil and guide wire for damage after shaping and before using.
- Use care when shaping the tip of this guide wire. Be sure the guide wire is wet before shaping to avoid damaging the hydrophilic coating.
- Verify which is the distal end before insertion and be sure to insert the flexible distal end (coiled end).
- Fasten the torque device to the guide wire firmly so that it may not loosen. Torque operation with a loose torque device may cause the coating to come off.
- When changing the attachment position of the torque device on the guide wire, loosen the fastening of the torque device before moving it.
- If any resistance is felt during insertion into a catheter used with the guide wire, do not use it.
- Do not manipulate the stopcock of the guiding catheter when the guide wire and/or other device is inserted in the guiding catheter fitted with a stopcock.
- Do not use a drug containing an organic solvent or an oil-based contrast medium, with this guide wire.
- Do not immerse the surface of the guide wire in alcohol, chlorhexidine gluconate solution, etc., or wipe it with gauze, absorbent cotton, etc. soaked with them. In addition, do not wipe the guide wire with dry gauze, absorbent cotton, etc.
- Do not repeatedly bend and stretch the same position of the guide wire, or continuously turn it in a curved vessel for a long period of time.
- Do not manipulate the guide wire with a torque device other than the supplied torque device. When attaching a torque device to the guide wire, or manipulating the guide wire using a torque device, carefully manipulate the guide wire so that the guide wire will not be damaged. Fasten the torque device with care not to excessively fasten it.
- To exert the lubricity of the guide wire's surface, inject heparinized saline into the guide wire holder tube and the lumen of the catheter to be used before using the guide wire. Confirm that the whole surface is wet, then take out the guide wire, and insert it into the catheter. To prevent eventual thrombus caused by the manipulation of the catheter and the guide wire, flush the catheter lumen with heparinized saline or other

suitable solution.

- Do not let the guide wire come into contact with liquids other than an intravascular contrast medium or heparinized saline.
- If bending is observed in the distal end of the guide wire during the procedure, do not manipulate the guide wire with the distal end being bent.
- Do not manipulate the guide wire while it is tightly secured with a Y-connector.
- If manipulation using two or more guide wires is performed, either guide wire may be tangled with the other guide wire or the other device inserted in the same catheter. Therefore, if it is necessary to perform manipulation using two or more guide wires, take extreme care.
- If the guide wire is bent or deformed due to an incidental load during use, stop using it.
- When manipulating this guide wire, check the movement of the product under high-resolution fluoroscopy. When moving this guide wire or the device used with the guide wire, always pay attention to their movement so that this guide wire or the device used with the guide wire does not interfere with each other or damage the vessel.
- If the guide wire meets resistance with the device used with the guide wire, do not apply an excessive force. If there is abnormal resistance, remove the whole system from the patient's body, and check the cause of the resistance.
- Remove blood or a contrast medium adhering to the surface of the guide wire using heparinized saline.
- Take preventive measures against infection after use. Discard this guide wire as medical waste.

Malfunctions and Adverse effects

During use of this guide wire, the following malfunctions and adverse effects may occur. If the malfunctions and adverse effects are serious, it may induce death or serious complication(s). Note, however, that malfunctions and adverse effects are not limited to these.

1. Malfunctions
 - Failure to cross a lesion
 - Breakage or bending of the guide wire
 - Damage such as separation
 - Difficulty in removal
 - Coming off of coating
2. Adverse effects (in random order)
 - Death
 - Infection
 - Vessel dissection
 - Damage to a vessel, including possible vessel perforation
 - Bleeding complications
 - Embolism
 - Thrombus
 - Infarction
 - Embolization or vascular damage caused by residue
 - Ischemia
 - Arrhythmia
 - Vasospasm
 - Vascular occlusion
 - Aneurysm rupture/perforation
 - Hemodynamic compromise
 - Aneurysm (false/dissecting)
 - Blood pressure reduction
 - Allergic reaction

How to Use

1. Inspection prior to use
 - 1-1. Before use, inspect carefully and confirm all devices and packages are undamaged.
2. Preparations for use
 - 2-1. Remove the holder tube containing the guide wire from the sterile pack.
 - 2-2. Before pulling the guide wire out of the holder tube, flush it with heparinized saline from the flush connector end into the holder tube to moisten the guide wire at least for 30 seconds. If it is difficult to pull the guide wire out of the holder tube, flush it again with heparinized saline.
 - 2-3. After pulling the guide wire out of the holder tube, inspect it to make sure that it is not damaged.
 - 2-4. Shape the distal end of the guide wire into the intended shape using the attached shaping device as needed. Inspect the guide wire for damage before and after shaping the distal end.
3. How to Use
 - 3-1. Before inserting the guide wire into an interventional device, wet it completely with heparinized saline.
 - 3-2. Insert the guide wire into an interventional devices using the accompanying inserter.
 - 3-3. Attach the accompanying torque device to the guide wire if necessary.
 - 3-4. Select the target vessel by carefully advancing or rotating the guide wire by fluoroscopy.
 - 3-5. After pulling the guide wire out of the body, remove blood and so on, and keep it wet.

Storage method

Do not keep the product in a bent and/or heavily-loaded condition. This product must be kept out of water. Store in a cool, dark and dry place.

Expiration date

The expiration date is indicated on the label of the device package.

Contents

1 set per package

Liability Disclaimer

By no means shall “ASAHI INTECC CO., LTD. and its affiliated companies” (hereinafter referred to as the “Company”) be liable for accidents, personal injuries, adverse effects due to any improper use of the product(s) or any other use inconsistent with these instructions. In no event shall the Company be liable for any damages either (i) arising out of storage of the product(s) after the shipment from the Company or (ii) due to selection of patients, surgery techniques, or any other medical activities by the medical institution that uses the product(s).

“ASAHI” is a trademark or a registered trademark of ASAHI INTECC CO.,LTD. in Japan and other countries.

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©2018 ASAHI INTECC CO., LTD.