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Edwards

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
Femoral Venous Cannulae
Pressure Drop vs. Flow Rate

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Manufacturer 
 Edwards Lifesciences LLC
 One Edwards Way
 Irvine, CA 92614 USA

949-250-7500
 800-424-3278

Edwards Lifesciences Services GmbH
 Edisonstr. 6
 85716 Unterschleißheim
 Germany

08/2018
 DCC-0055472 C
 © Copyright 2018, Edwards Lifesciences LLC
 All rights reserved.

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Edwards

FEMORAL ACCESS VENOUS AND ARTERIAL CANNULA

Instructions for Use

Directory

English (EN)

English (EN)

FEMORAL ACCESS VENOUS AND ARTERIAL CANNULA

Instructions for Use

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

The Edwards Lifesciences Femtrak, VEM, and Fem-Flex II femoral venous cannulae and Fem-Flex and Fem-Flex II femoral arterial cannulae are wire-reinforced, thin-wall cannulae. The clear proximal section is not reinforced to allow clamping and terminates in a barbed connector. Each cannula is supplied with a dilator and/or introducer that will pass over a 0.038" (0.96 mm) guidewire. Barium stripes on the tubing aid in cannula placement. Some codes include an additional non-lumen dilator. Some product codes may also include a suture ring. The Fem-Flex II (FEMIIA) cannulae include two barium stripes on the tubing and an additional non-lumen dilator to aid in placement.

Each Edwards Lifesciences device is packaged sterile and non-pyrogenic in a sealed, peel-type pouch.

INDICATIONS FOR USE

Edwards Lifesciences Femoral Access Cannulae are intended for use in situations which require rapid femoral venous and arterial access for short-term (≤ 6 hours) cardiopulmonary bypass. Vessel access (venous or arterial) is left to the discretion of the physician.

Femoral access cannulae may be used in pediatric populations or adult populations based on flow rate requirements and individual patient anatomy. Please consult labeling to determine pressure drop related to flow rates.

CONTRAINDICATIONS FOR USE

Femoral Access cannulae are contraindicated for long-term use (> 6 hours), including Extra Corporeal Membrane Oxygenation (ECMO) procedures.

Devices are contraindicated for placement in vasculature outside of the femoral artery or vein.

Single tube venous drainage from the right atrium should not be used in procedures where air from the procedure site may enter the cannula.

This device is not intended for use other than as indicated and should not be used when any physical impairment would contraindicate its use.

WARNINGS AND PRECAUTIONS

Supplied sterile and non-pyrogenic in undamaged package. Do not use if device shows signs of damage (i.e., cuts, kinks, crushed areas, leakage), or if package is damaged or open as this may indicate compromised sterility and/or product damage.

Use fluoroscopy and/or ultrasound imaging during placement to ensure proper positioning and product performance.

Failure to properly advance the introducer/cannula or introducer/sheath over a previously advanced guidewire may result in vessel perforation and/or dissection.

FOR SINGLE USE ONLY. This device is designed, intended, and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE THIS DEVICE. There are no data to support the sterility, non-pyrogenicity, and functionality of the device after reprocessing. Such action could lead to illness or an adverse event, as the device may not function as originally intended.

Products are known to contain phthalates, which may be found in device materials containing plasticizers such as DEHP and BBP. High exposure to such phthalates during medical treatments in children and pregnant or nursing women may raise a concern. A review of available data and literature supports the conclusion that the benefits outweigh the overall residual risk.

Proper surgical procedures and techniques are the responsibility of the medical profession. Described procedures are provided for informational purposes only. Each physician must determine the appropriate use of this device for each patient based on medical training, experience, the type of procedure employed, and the benefits and risks associated with device use.

This device is intended for short-term use only (≤ 6 hours). There are no data to support the function and performance of the device beyond six hours of use.

Use the largest cannula possible, consistent with good surgical technique, to ensure adequate venous drainage.

Surgical entrance or introduction of air into the right heart may result in the partial or total loss of venous drainage to the extracorporeal circulation.

If increased resistance is felt at any time upon insertion or removal of the guidewire, dilators, introducer, or cannula, stop and investigate the cause before continuing. Inability to easily advance or remove these devices may indicate vascular disease or injury. Closely examine the device position within the vessel using fluoroscopy and/or ultrasound imaging prior to proceeding.

Avoid over-insertion of the cannula to reduce risk of cannula tip impingement against the vessel walls which may lead to excessive line pressures and vessel damage.

Avoid over-insertion of the cannula to reduce risk of cannula tip impingement against the vessel walls which may lead to decreased drainage and vessel damage.

Over-insertion of the cannula may result in proximal holes entering the superior vena cava. This may lead to diminished venous return. Cannula markings and/or palpation is recommended to ensure proper cannula placement.

Line pressure exceeding acceptable clinical limits may result from incorrect cannula tip position and/or restricted tip patency. Confirm proper placement to minimize risk of injury.

Securely tie-band the connector to cannula junction prior to initiating bypass to protect against tubing slippage.

Insert femoral cannula using standard surgical techniques, including Seldinger technique and cut-down under direct visualization, to ensure proper placement.

Ensure proper levels of anticoagulant therapy are maintained prior to insertion of the cannula and throughout cardiopulmonary bypass, to reduce the risk of complications due to thrombus formation on or within the cannula, and in the blood stream.

When provided, ensure caps are securely seated to reduce the risk of leakage.

Ensure that the vessel used is of adequate size to allow sufficient perfusion distal to the cannula after insertion.

Reinsert introducer only if tip is visible to avoid potential vessel damage.

Wire-reinforced cannulae should be clamped in the non-reinforced section located at the connector end since clamping of the reinforced section may produce permanent cannula deformation, thereby impeding flow through the cannula and risking puncture or tearing of the cannula.

The leg used for femoral drainage/perfusion may receive inadequate collateral circulation resulting in severe ischemia. CLOSELY MONITOR THE INVOLVED EXTREMITY FOR ANY SIGN OF DEVELOPING ISCHEMIA.

For cut-down procedures, directly visualize vessel during insertion and removal of the cannula to minimize vessel damage.

Monitor vessel patency and hemostasis after removal of cannula to avoid possible ischemia and/or blood loss.

When using small diameter cannulae, do not exceed the maximum pressure limit while increasing flow. A rapid pressure rise may result.

Prior to initiating cardiopulmonary bypass, completely expel air from the system. Failure to eliminate air from the system could result in an air embolus.

Dispose of used product in accordance with established hospital protocols for biohazards to minimize risk of exposure.

WARNINGS AND PRECAUTIONS SPECIFIC TO FEMTRAK

Not recommended for use with coated guidewires as they may not be compatible with this product.

DIRECTIONS FOR USE

WARNING: When using small diameter cannulae, do not exceed the maximum pressure limit while increasing flow. A rapid pressure rise may result.

1. Use fluoroscopy and/or ultrasound imaging during advancement, positioning, and placement of the cannula.
2. Cannulate for femoral bypass using standard surgical techniques, including Seldinger technique and cut-down under direct visualization to ensure proper placement.
3. If insertion components are not included, Edwards' PKV or PIKA kit is recommended.

Note: A moveable suture ring is available with some product codes.

CAUTION: Ensure that the vessel used is of adequate size to allow sufficient perfusion distal to the cannula after insertion.

4. Insert cap/introducer assembly into cannula. Seat cap securely on connector. Following placement and removal of introducer, the cannula may be clamped in the area between the wire reinforcement and the connector. Remove cap from the end of the connector before connecting to the extracorporeal circuit. Graphs of pressure-flow relationships are provided at the end of these instructions.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Edwards, Edwards Lifesciences, the stylized E logo, Fem-Flex, Fem-Flex II, Femtrak, and VEM are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

Made in the U.S.A.


CHART LEGEND

	Femoral Arterial Cannulae Pressure Drop vs. Flow Rate	Pressure Drop (mmHg)	Flow Rate (L/min), H ₂ O at Room Temperature
--	---	----------------------	---


CHART LEGEND

	Femoral Venous Cannulae Pressure Drop vs. Flow Rate	Pressure Drop (mmHg)	Flow Rate (L/min), H ₂ O at Room Temperature
--	---	----------------------	---

SYMBOL LEGEND

EN	Single use	Sterilized Using Ethylene Oxide	Non-pyrogenic	Do not use if package is damaged	Lot Number	Use By	Authorised Representative in the European Community	Date of Manufacture

					Rx Only	
EN	Consult instructions for use	Caution	Quantity	Contains phthalates	Catalogue Number	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Note: Not all symbols may be included in the labeling of this product.

Edwards

Manufacturer 
Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

EC REP

Edwards Lifesciences Services GmbH
Edisonstr. 6
85716 Unterschleissheim
Germany

08/2018
DOC-0055472 C
© Copyright 2018, Edwards Lifesciences LLC
All rights reserved.

949-250-2500
800-424-3278

DATA MATRIX
BARCODE FPO