



ASAHI CHIKAI^誓 X 010

Neurovascular GuideWire


ASAHI CHIKAI^誓 X 010

Neurovascular Guide Wire

~Gold standard in distal vessel cases~

By enhancing flexibility, reduction of resistance from the vessel wall is enabled and superior torque response is achieved.

Structure and Ordering Information



Product	Catalog No.	O.D.	Coating	Usable Length	Coil Length	Radiopaque Length	Tip Shape
ASAHI CHIKAI X 010	AIN-CKX-10-200-R	0.26 _{mm} (0.010 _{inch})	Hydrophilic coating (SLIP-COAT®) 180 _{cm}	200 _{cm}	16 _{cm}	16 _{cm}	Re-shapeable round curve

* Re-shapeable angle 25°

CE 0344

© 2019 ASAHI INTECC CO., LTD.
ASAHI, and *CHIKAI* are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries.
Coated with SLIP-COAT® coating.

Manufactured by:

ASAHI INTECC

ASAHI INTECC CO.,LTD.

3-100 Akatsuki-cho, Seto, Aichi 489-0071 Japan
TEL:+81-561-48-5551 FAX:+81-561-48-5552

www.asahi-intecc.com

Distributed by: