

# SAPPHIRE<sup>®</sup> II PRO

Coronary Dilatation Catheter

## Thông tin đặt hàng


### EU Manufacturing Catalogue Numbers

Đường kính bóng (mm)	Chiều dài bóng (mm)				
	10	12	15	20	30
1.75	217-103-5U	-	217-153-5U	217-203-5U	-
2.0	220-103-5U	220-123-5U	220-153-5U	220-203-5U	-
2.25	222-103-5U	-	222-153-5U	222-203-5U	-
2.5	225-103-5U	225-123-5U	225-153-5U	225-203-5U	225-303-5U
2.75	227-103-5U	-	227-153-5U	227-203-5U	-
3.0	230-103-5U	230-123-5U	230-153-5U	230-203-5U	230-303-5U
3.25	232-103-5U	-	232-153-5U	232-203-5U	-
3.5	235-103-5U	-	235-153-5U	235-203-5U	235-303-5U
4.0	240-103-5U	-	240-153-5U	240-203-5U	-

## Bảng thông số độ giãn nở

Áp lực bơm (atm)	Đường kính bóng (mm)								
	1.75	2.0	2.25	2.5	2.75	3.0	3.25	3.5	4.0
2	1.66	1.87	2.11	2.36	2.59	2.82	3.06	3.29	3.77
4	1.70	1.94	2.18	2.43	2.67	2.91	3.16	3.39	3.89
6 NOM*	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	4.00
8	1.80	2.06	2.32	2.57	2.83	3.09	3.34	3.61	4.11
10	1.84	2.13	2.39	2.64	2.91	3.18	3.44	3.71	4.23
12	1.89	2.19	2.45	2.71	2.98	3.27	3.53	3.82	4.34
14 RBP**	1.94	2.25	2.52	2.78	3.06	3.36	3.62	3.92	4.45
16	1.98	2.32	2.59	2.85	3.14	3.45	3.72	4.03	4.57
18	2.03	2.38	2.66	2.92	3.22	3.54	3.81	4.14	4.68

\* Nominal Pressure. The nominal in-vitro device specifications do not take into account any lesion resistance.  
 \*\* Rated Burst Pressure. Do not exceed RBP.



For more information please visit our website at [www.OrbusNeich.com](http://www.OrbusNeich.com) or contact us:

Corporate Headquarters  
 Unit 303 & 305, 3/F, Building 20E  
 Hong Kong Science Park  
 Shatin, N.T., Hong Kong, China  
 Phone +852.2802.2288  
 Fax +852.2507.3532

EMEA Regional Headquarters  
 Drs. W. van Royenstraat 5  
 3871 AN Hoevelaken  
 The Netherlands  
 Phone +31.33.254.1150 Toll Free Phone 00800.0254.1150\*

\* Only for Belgium, Denmark, France, Germany, Ireland, Netherlands, Norway, Sweden and UK

©2020 OrbusNeich Medical Company Limited or its affiliates. All rights reserved.  
 The Sapphire II PRO product is now approved in the USA and distributed with USA specific catalog numbers.

G-70-0531 Rev03

Sapphire<sup>®</sup> and OrbusNeich<sup>®</sup> are registered trademarks of OrbusNeich Medical Company Limited or its affiliates.  
 All data and photos on file. Illustrations are not to be considered as engineering drawings or photos.

 **OrbusNeich<sup>®</sup>**  
 Pioneers in life-changing technologies

# SAPPHIRE<sup>®</sup> II PRO

Coronary Dilatation Catheter



Tiêu chuẩn mới cho  
**Bóng nong mạch vành**

 **OrbusNeich<sup>®</sup>**  
 Pioneers in life-changing technologies


## Tự tin vượt tổn thương Mọi thời điểm

Bóng Sapphire II PRO được thiết kế đặc biệt nhằm vượt mọi vết thương khó nhất và những tổn thương nghèo nàn, đáp ứng mọi kích thước của tổn thương




## Hệ thống dẫn với công nghệ P-TECH cho lực đẩy tối ưu và chống xoắn

Hệ thống dẫn với công nghệ P-TECH được thiết kế độc nhất với ống thông tiết diện nhỏ có khả năng chống xoắn và khả năng vận chuyển xuất sắc




## Công nghệ đầu tip Sub-zero giúp vượt qua tổn thương trơn mượt

Đầu tip thuôn nhọn cân bằng tốt với công nghệ Sub-Zero có tiết diện cực nhỏ, đảm bảo khả năng đâm xuyên những tổn thương nhỏ nhất một cách dễ dàng nhất



Emerge

Sapphire II PRO



Lacrosse LAXA


Sapphire II PRO



## Bóng XR được thiết kế cho những trường hợp bất khả thi

Bóng XR độc quyền với khả năng đâm xuyên và thu hồi tốt nhất mà không ảnh hưởng đến độ bền và độ mạnh

So sánh đường kính thân bóng với bóng có đường kính 3.0 mm x 15 mm



\*Data on file, scale depictions